

pliance at a glance

■ ftm analyser

number of sensors (max.)	12544
offset and gain	programmable for each sensor
scanning speed	10000 sensors / second
synchronisation in and out	TTL, selectable
analog input channels (max)	21
wireless interface	available
power supply	12 V DC or 110 to 240 V AC

■ ftm sensor mat (standard)

technology	capacitive matrix sensor
number of sensors	256 ~ 1344
sensor resolution	2.0, 2.5, 6.0 cm ² / sensors
pressure range	1 ~ 60 kPa
hysteresis	< 3 % FSO
temperature coefficient	0.05 kPa x °C
elasticity	approx. 4 %
min.bending radius	40 mm

novel_{gmbh} (Germany) • Ismaninger Str. 51 • D-81675 München • Tel: (+49) 89-417767-0 • Fax: (+49) 89-417767-99 • e-mail: novel@novel.de
 novel_{electronics} inc. (USA) • 964 Grand Avenue • St.Paul, MN 55105 • Tel: (+1) 651-221 0505 • Fax: (+1) 651-221 0404 • e-mail: novelinc@novel.de
 novel_{gmbh} (Great Britain) • Tel/Fax: (+44) 20-8659 0959 • e-mail: noveluk@novel.de

www.novel.de

emed[®], pedar[®], pedograph[®], pliance[®] and the Logo (colored foot) are registered trademarks of novel _{gmbh} © 2002

Information subject to change without notice



pliance_ftm_11/06/2002_mindae

pliance system

pliance systems offer the state of the art technology for pressure distribution measurement between soft and curved surfaces.

The system consists of a flexible and elastic measuring mat, a multi-channel analyser, a calibration device and a software package for PC's.

The elastic measuring mats are available in various sizes, sensor configurations and pressure ranges.

The pliance-ftm analysers vary from small portable 16x16 channel units to large 112x112 channel units with a wide range of options, such as master-slave synchronisation of several systems, dynamic amplification control, synchronisation of video systems and analog inputs for accelerometers.

The patented **trublu**® calibration device can be used at any time to verify the quality of the measuring results. It works with homogeneous air pressure on all sensors through increasing steps of pressure. Thus an individual calibration curve for each sensor is calculated and used during data acquisition.

The software operates as a Windows application and contains many useful methods of data collection and scientific analysis of dynamic pressure distribution. It also allows continuous data storage in online mode and data handling with a configurable SQL database. The expert can design the parameter configuration to specific needs and exchange data with colleagues via HTML protocol. As in all physical measuring

systems the most important part is the sensor technology. pliance works with capacitive transducers in a matrix configuration. The elasticity of the sensor mats permits perfect conformability to 3-dimensional deformations.


The pressure transducing elements contain high-tech elastomers manufactured by novel. Restoring force, range of force, threshold, hysteresis, temperature effect, frequency response and other characteristics are determined during the manufacturing process. This makes it possible to adapt the sensor characteristic to different measuring needs.




trublu® calibration device

novel developed an analyser technology that allows not only individual calibration curves for each sensor, but also individual dynamic amplification control and crosstalk suppression, resulting in very accurate and reproducible pressure values.


pliance software




Icon selection



Online display for record and play



3 - D and x, y cross - section display




Isobar display with horizontal cross sections

novel database pro 14 - C:\nov1026\novfile\database\VBPRO14\demo_carseat.gdb

ID number	seat type	car model
0	novel Demo	
1	Seat 1	Porsche
2	Seat 2	BMW
3	Seat 3	Honda

name	sex	body mass	height
CF	male	65	179
DJ	female	55	170

testnumber	date	exper
1	15.10.2002	static
2	15.10.2002	apply



novel report

